



FRIEDRICH-SCHILLER- UNIVERSITÄT JENA NMR Platform

- **How to download the free (academic)
TopSpin program -**



Home - About us - Register

Register for Bruker users account

Welcome to Bruker quick registration portal. To download Bruker software, you are required to have an account. If you don't have an account with us yet, you can start now by filling out the form below.

Please also note that the personal information that you enter in the form below may be screened against relevant sanction lists to ensure compliance with applicable export control laws.

First Name: *	<input type="text"/>	I am interested in the following newsletters: <input type="checkbox"/> Magnetic Resonance eNews & The Resonance Newscast <input type="checkbox"/> Molecular Vibrational Spectroscopy <input type="checkbox"/> Mass Spectrometry <input type="checkbox"/> FIRST Newsletter (X-ray Diffraction and Elemental Analysis) <input type="checkbox"/> Nanovations (Atomic Force Microscopy) <input type="checkbox"/> Surface Insights (Surface Profiling and Mechanical Testing) <input type="checkbox"/> Preclinical Imaging e-news <input type="checkbox"/> BioViews (Life Sciences Microscopy) <input type="checkbox"/> CBRNE Detection e-news <input type="checkbox"/> Microbiology & Diagnostics
Last Name: *	<input type="text"/>	
Country or Area: *	Germany	
Email Address: *	test@uni-jena.de	
Create Password: *	<input type="password"/>	
Confirm Password: *	<input type="password"/>	
Password must be at least 6 characters long		
<input type="checkbox"/> I'm not a robot		
<input type="button" value="Register"/>		

Privacy Information

Please view our [Privacy Notice](#).



For using the **free academic version** of TopSpin (with a duration of up to 15 years!!) create a user account on the Bruker webpage, using your email adress of the university and a new password.

(<https://www.bruker.com/about-us/register.html>)



Login to Bruker

Home - Login

Benutzeranmeldung

Enter your username and password here in order to log in on the website:

Email Address

Password

[Reset password](#) [Register for an account](#)
[Request for new authentication link](#)




Log in ...



Home - Service - Support und Upgrades - Software-Downloads - NMR - NMR TopSpin License for Academia

NMR TopSpin License Processing for Academia

In order to provide students, researchers and teachers with unlimited access to the best tools for off-line NMR processing, Bruker is making their market leading NMR processing software TopSpin available free of charge for all academic users. In addition to Bruker NMR data, the software is able to process data acquired on Agilent, Varian and Jeol systems, as well as read any data provided in JCAMP format. We are happy to offer you a free Academia license for **TopSpin 4 and TopSpin 3.6** ("Academia license").



Disclaimer


Access to Academia licenses is restricted to academic customers. By using this web page for license requests you confirm that you are employed by or studying at an academic- or governmental institution and that you use the software exclusively in this context. Your request and the license generation will be logged. If you do not agree with this, please do not continue.

In the next step a personal ticket code will be generated for you. With this ticket code you can activate the Academia license on a compatible computer of your choice and use the software:

The use of Academia licenses for commercial organizations or companies is prohibited. To order licenses for a commercial organization, please follow this link.



... and choose your operating system.



NMR Software Downloads for Windows

Kontakt

Home - Service - Support und Upgrades - Software-Downloads - NMR - PC - PC TopSpin

TopSpin - Windows PC

Current versions:

- TopSpin 4.0.6 (for AVANCE NEO and processing in general)
- TopSpin 3.6.1 (for AVANCE III, III HD and processing in general)
- TopSpin 3.2 pl7 (for AVANCE II, III, III HD, Fourier 300 and processing in general)
- TopSpin 3.1 pl7 (for AVANCE II, III, Fourier 300 and processing in general)
- TopSpin 2.1 pl8 (for AVANCE I)
- TopSpin 1.3 pl10 (for AVANCE DRX, DPX, DMX, and DSX)

For detailed information refer to the corresponding [release letter and software support matrix](#).

(For information on acquiring the current version of TopSpin on DVD, contact your nearest [Bruker BioSpin office](#)).

To work with TopSpin a license is required. To download the current TopSpin version, click on the link and follow the instructions provided below. Please store the file on your hard disk and start the installation when the download is complete.

Current Release Download

 [TopSpin 4.0.6 & CMC-assist 2.12 \(1663.9 MiB\)](#)

(Created: 2018-11-21)
TopSpin for AVANCE NEO and processing in general.
This link allows for the direct installation of a complete acquisition/processing version for Windows 7 and Windows 10.
MD5 checksum:
d59c742cfbd73dccc1b767bf12e5287


ISO - Downloads

These images contain everything that you would also see on a DVD when you order it from Bruker. They contain all programs for all supported operating systems and also all additional software like AMIX, AURELIA and CMC-assist. For programs that are supported on the MAC, the MAC version is also contained in this image.

 [TopSpin 4.0.6 \(rel40.iso\) \(6175.1 MiB\)](#)



Download the current release version of TopSpin.



After the installation (see 'Installationanleitung TopSpin' on our NMR Platform webpage) request your personal license ticket.

Free TopSpin Processing for Academia

Kontakt

Home - Service - Support und Upgrades - Software-Downloads - NMR - NMR TopSpin License for Academia

NMR TopSpin License Processing for Academia

In order to provide students, researchers and teachers with unlimited access to the best tools for off-line NMR processing, Bruker is making their market leading NMR processing software TopSpin available free of charge for all academic users. In addition to Bruker NMR data, the software is able to process data acquired on Agilent, Varian and Jeol systems, as well as read any data provided in JCAMP format. We are happy to offer you a free Academia license for **TopSpin 4 and TopSpin 3.6** ("Academia license").

Free License For Academia

With the following code you can activate your free academia license:

ABCDE-12345-FGHIJ-678EX-AMPLE

Just write it down (or copy it to the clipboard), start TopSpin, choose the option "Activate this installation now" and follow the instructions. If you have not installed TopSpin already you can download it [here](#).

Your license is valid for 15 years. You can use it on a compatible computer of your choice and you can move it to another system if needed. See the section "Rehosting of license" in the license management manual for details.

[Download for Windows](#) [Download for Linux](#) [Download for Mac OS](#)

Disclaimer

Access to Academia licenses is restricted to academic customers. By using this web page for license requests you confirm that you are employed by or studying at an academic- or governmental institution and that you use the software exclusively in this context. Your request and the license generation will be logged. If you do not agree with this, please do not continue.

In the next step a personal ticket code will be generated for you. With this ticket code you can activate the Academia license on a compatible computer of

This website uses cookies for analytics and functionality purposes. To change your cookie settings or find out more, [click here](#). If you continue browsing our website, you accept these cookies.

OK

Use your personal code for activating the license within the installed software and follow the instructions.

Contact me for further question (f.pielenz@uni-jena.de).
Good luck!