

TopSpin Download Instructions

- > Deinstallation of older versions of TopSpin
- > Registration at the Bruker website (<https://www.bruker.com/en/legal/register.html>) with your FSU-mail address
- > Validation of said mail address
- > Log in on the Bruker website: <https://www.bruker.com/protected/en/services/software-downloads/nmr/nmr-topspin-license-for-academia.html>
- > „Request Academia License" and copy the license
- > Download TopSpin for Windows (<https://www.bruker.com/protected/en/services/software-downloads/nmr/pc/pc-topspin.html>) --> **TopSpin 4.1.4 & CMC-assist 2.21** (2088 MB) (for Win10)
- > Open the installation file and follow the instructions on screen
- > Once asked, paste the copied academic license from earlier

If any problems arise during the installation of the program, please feel free to contact the NMR Platform at **03641 9 48258** or nmr-service@uni-jena.de

